



KANOMAX

The Ultimate Measurements



Handheld Condensation Particle Counter

Model 3800



Applications

- Indoor Air Quality Investigation
- Occupational Health Control
- Gasket Leakage Testing
- Engine Exhaust Testing
- Filter Testing
- Aerosol Research

Features

- Particle size range: 0.015 to 1 μ m
- Measuring range: 0 to 100,000 particles/cm³
- Sampling airflow: 0.7L/min
- Data logging up to 100,000 measurements
- Measuring Software included (for Windows)
- USB interface for data transfer
- Power supply: 6 x AA batteries or AC adapter
- Weight: 3.3lbs (1.5kg)

Benefits

- 0.015 micron (15nm) minimum sensitivity
- Handy and Light weight
- Data storage and logging capability
- Simple to download the data to your PC via USB

Optimal Screening Tool for Nano Particles

■ Applications



Checking worker exposure to airborne contaminants



Leakage Testing



Engine Exhaust Testing

